F LCS Perox

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GB LCS Perox

LCS Perox is a very modern, simple and effective system for daily contact lens care.

By combining the disinfecting power of hydrogen peroxide enhanced by the simultaneous action of a carefully-formulated tablet, the cleaning effect is greatly increased (and selectively targets protein deposits) and the peroxide is neutralised slowly and completely to give a suitable saline solution:

- Neutral pH
- Isotonic with tear
- No preservatives. Accordingly, lenses are kept perfectly clean, disinfected, hydrated and ready for use maintaining their original parameters.

The **LCS Perox** system contains:

- Bottle of solution.
- Blister pack of tablets.
- Lens case with elastic valve.

INSTRUCTIONS FOR USE:

1. Fill the lens case up to the mark with **LCS Perox** Solution.
2. Place the lenses in their respective lens holders.
3. Add one tablet in the lens case and screw the lid on.
4. Shake the lens case for a few seconds.
5. Never use the solution without the neutralizing tablet.
6. Wait at least 2 hours for the cleaning, disinfection and neutralisation process to conclude.
7. The lenses can be used after 2 hours or can be left in the lens case overnight.
7. The neutralising tablet contains chlorophyll. When the process is finished, the solution where the lenses are should be light green, which indicates that the tablet has worked.
8. Shake the lens case for a few seconds.
9. Before handling the lenses to fit them, wash your hands and rinse well to eliminate all traces of soap.
10. Empty the solution out of the lens case.
Keep the empty lens case closed when not in use.

Last revision date: 2009/11

PRECAUTIONS:

- Take all hygiene precautions before handling your contact lenses (for example, wash your hands).
- Avoid prolonged contact of the initial solution with your skin.
- If necessary, rinse with tap water.
- Do not use the product if the safety seal is missing or broken.
- Do not ingest the solution and keep out of children's reach.
- Once opened, the bottle should not be used after 60 days.
- Maintain the bottle well closed after use.
- Once the cleaning and disinfection process is completed, the lenses are in a sterile solution WHICH CONTAINS NO PRESERVATIVES.
- If the lenses are not used on a daily basis, the green colour of the solution will fade after a few days. Because it does not contain preservatives, if the solution is colourless, this is a sign that you should disinfect the lenses again before wearing them.
- If you do not use your lenses every day, it is advisable to replace the cleaning solution in the lens case with a solution containing preservatives, for example **LCS Saline**, to avoid possible contamination after the disinfection process has concluded.
- **LCS Saline**, may also be used for rinsing your lenses (if they fall, or are temporarily removed, etc), instead of carrying out a full disinfection with **LCS Perox**.
- If, exceptionally, the lens wearer's physiology leads to considerable protein deposits that **LCS Perox** is unable to remove, seek the advice of your contact lens specialist regarding the advisability of additional protein removal (**LCS Enzym** is suitable for this purpose).
- Discard the solution in the lens case after each use.
- Never use tap water to rinse your lenses or lens case.
- Keep at room temperature.
- In case of irritation, consult your contact lens specialist.
- Before using any ophthalmic medicine, inform your doctor you are a contact lens user.
- Do not use after the expiration date.
- Do not use if you are allergic to any of the components.
- Do not put unneutralised disinfectant solution directly in your eye.

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