



EasySoft[®] Confort

Hydrogel Soft Silicone lens with
biannual renewal for each ametropia

Indications

EasySoft Confort is a biannual soft Silicone Hydrogel lens for each ametropia :

- **Myopic and hyperopia**
EasySoft Confort : for all spherical ametropias
- **Astigmatism**
EasySoft Confort Toric : Anterior toric

Technical information

Matériau	Silicone hydrogel Filcon v3 Other material on request)
Radius R_0	From 7,80 to 10,50 in 0,10 mm steps
Diameter \varnothing_T	From 12,00 to 16,00 mm in 0,10 mm steps
Power P	From + 30,00 to - 30,00 δ in 0,25 δ steps Cylinder : From - 0,50 to - 6,00 δ in 0,25 δ steps Axis : all axis par 5°
Packaging	Single bottle

Wear and care

Wear	Daily
Care	Oté Clean and multifunction for soft lenses

Adaptation

EasySoft Confort

$\varnothing_T = 14,00$ mm

$R_0 = 8,60$ mm for $7,60 \leq K \leq 8,20$ mm

P= sphere VL brought back to the corneal top

EasySoft Confort Toric

$\varnothing_T = 14,50$ mm

$R_0 = 8,60$ mm for $7,60 \leq K_m \leq 8,20$ mm

P = sphere VL brought back to the corneal top

Tél. +33 231 86 66 56
info@laboratoire-lcs.com

